

	<p>BagMixer® 400 CC® 400 ml lab blender</p> <p>Ref. 024 230</p> <p>Powerful and Click and Clean®. Avoiding all risks of cross-contamination, BagMixer® 400 CC® is an easy-to-use and powerful lab homogenizer. Adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, BagMixer® 400 CC®, equipped with its window door and patented adjustable paddles, ensures fast and perfect homogenization for lab analyses. Cleaning is easy with Click and Clean® removable paddles system.</p>
---	---

TECHNICAL SPECIFICATIONS

- ▶ All stainless steel
- ▶ Window door
- ▶ Multifunction digital screen
- ▶ Silent brushless motor
- ▶ Side-by-side paddle stop
- ▶ Useful volume: 50 – 400 ml
- ▶ Variable speed (4, 6, 8, 10 strokes/sec.)
- ▶ Variable time (1 sec. – 1 h or ∞/count down)
- ▶ Variable power, adjustable paddles (patented) with on-screen indication
- ▶ Q-Tight® closing: no more upward leaks
- ▶ Click & Clean® removable paddles system (patented)
- ▶ Fully opening door (270°)
- ▶ Security drip tray

- ▶ Power: ~100-240V / 50-60Hz
- ▶ Machine: size (L x W x H): 42 x 26 x 26.5 cm, weight: 16 kg
- ▶ Box: size (L x W x H): 48 x 33.5 x 33.5 cm, weight: 17 kg

- ▶ 3 year guarantee (after registration)
- ▶ Window door and shock absorbers lifetime guarantee
- ▶ Lifetime service
- ▶ In compliance with: ISO 7218 standard
- ▶ Designed and made in France
- ▶ Manufactured under: ISO 9001 V2000, CE, WEEE 2002/96/EC, RoHS 2002/95/EC standards

DELIVERED WITH

- ▶ Security drip tray
- ▶ Power cord
- ▶ User manual, IQ/OQ/PQ guidelines
- ▶ Quick user guide
- ▶ Free samples of sterile bags: BagPage® Plus 400, BagFilter® 400 P, BagFilter® 400 S

OPTIONS / ACCESSORIES

- ▶ Ref. 211 040 BagOpen® 400, bag opener stand for holding bags
- ▶ Ref. 221 040 BagRack® 400, storage racks for 10 bags
- ▶ Ref. 231 040 BagClip® 400, closing clips for lab blender bags 400 ml
- ▶ Ref. 251 091 BagPipet®, automatic pipet, fixed volume 0.1/0.9/1 ml, ejector for BagTips®

CONSUMABLES

- ▶ Ref. 111 425 BagFilter® 400, sterile side-filter bag
- ▶ Ref. 122 025 BagPage® 400, sterile full-surface filter bag
- ▶ Ref. 132 025 BagLight® 400, sterile plain bag
- ▶ Ref. 145 040 RollBag® 1300, sterile sampling bag
- ▶ Ref. 252 019 BagTips® REGULAR, 19 cm sterile tips for BagPipet® or other pipet
- ▶ Ref. 252 024 BagTips® JUMBO, 24 cm sterile tips for BagPipet® or other pipet