

Corning Incorporated
Life Sciences

Registered ISO 9001

Catalog Number PAS3A-01
PDS reference PDS_Cbor_23926 Version _ 2
Product description Corning® Gosselin™ Transfer Pipette, 3 mL, Graduated, PE, Sterile, 1/Bag, 200/Box, 2000/Case

Sanitary Assurance Level

This product has been treated by ethylene oxide.

Expiration Statement

This product is designed to have an expiration date of 2 years from the date of manufacture regarding sterility.

Performance testing (AQL)

Inspection based on a statistically significant number of samples towards to internal or agreed specifications that determine Assurance Quality Level.

Primary body composition

PE

Secondary body composition

HDPE

Closure

N/a

Special feature

PE

Packaging

1/Bag, 200/Box



Hugues Vanderschelden
Plant Quality Manager, Borre
vandersch2@corning.com